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**INSTRUCTION FOR USE
ONCO-GRIP® AND HUBER NEEDLE WITH EXTENSION LINE**

CE marking: 2000

Products references: GNXXXX - LRXXXX

Product sterilized by ethylene oxide

Sterile product as long as the packaging hasn't been opened damaged or brake. Single use product

WARNING : Don't resterilize

WARNING : Store in a cool and dry place with ambient temperature. Avoid extreme temperatures and humidity. Avoid rapid and repeated changes of temperature during storage.

Conserve in the original packaging.

WARNING : Read carefully instructions before use.

WARNING : Check the integrity of the sterility protector primary packaging. Never use a device which primary packaging is damaged.

WARNING : Shelf life – Use before expiration date indicated on the package.

WARNING : NOT safety medical device

INDICATIONS

The ONCO-GRIP® and the Huber needle with extension set are composed of a Huber needle connected to a 18 cm length extension cable ended by a Luer Lock hub. The extension cable contains a closing clamp.

Its movable wings permit to have a good handing of the needle and thus precision during insertion. Its weak thickness and "blunt tray" allow a simple adhesive dressing use for a better comfort of the patient. The risks of needle tear out are much reduced.

They are intended to be used with access ports for:

- injection or liquid perfusion

- parenteral nutrition

The cannula shape allows port access limiting to the maximal the septum damages and the patient's pain.

The integration of an extension cable with Luer Lock hub to the device bring advantages face to face a simple Huber needle:

- easy use (the hub is much more accessible)

- minimize the infection risk eliminating a connection and insuring a perfect watertightness during syringe change (clamp).

The ONCO-GRIP® has movable wings allowing a good handing of the needle and thus precision during insertion. Its weak thickness and "blunt tray" allow a simple adhesive dressing use for a better comfort of the patient. The risks of needle tear out are much reduced.

USING PRECAUTION :

This device is not made to be used more than 3 to 4 days.

It is advised to change it daily or more often. In case of replacement by an identical device, the total period of use should not exceed 29 days

WARNING : Only use with devices with continuous flow or under direct control of the operator (syringe, ...).

Only the specifically identified references are compatible with high-pressure injectors. Do not use references in this manual in combination with a high pressure injector. High pressure or use of powerful injection pumps can cause leakage or damage to the system.

WARNING : Huber needles should only be used by trained personnel.

Use preferably 22G (0,7 mm) diameter or least for all the applications which don't require a big flow (injection, medicinal perfusion, heparinisation, rinsing...) in order to protect the septum of the access port. For the big flow applications (parenteral nutrition, administration of blood or of blood derivatives, it is preferable to use Huber needles with 20G (0,9 mm) diameter or superior according to the case.

CAUTION: Use only 10mL syringe minimum to prevent excess pressure in the system.

CAUTION: Never use an infected area.

CAUTION: Carefully tighten all connections before and after use. Failure to attach caps can result in an embolism or bleeding.

NEEDLE INTRODUCTION :

WARNING: The biggest careful has to be brought on the manipulation's ASEPSIS. Check that the introduction area doesn't present any infection sign (red blotches, local inflammation or pain at palpation).

WARNING: Implantable port puncture: slightly move up from the last puncture area in order to avoid elulice or fistula between the septum and the patient skin.

WARNING: Check that the length of the needle is correct with respect to the depth of the port body. Select a needle size adapted to the depth of the port body and to the thickness of the tissues. Too long needle can damage the access port; a needle too short may not fully go through the septum and lead to fluid injection from the port and / or blocking of the needle.

• Check that the needle is correctly placed in the port prior to injection. If there is a doubt about the placement of the needle, make an X-ray test to confirm that the positioning is consistent with the approved protocol.

WARNING: Do not use with a high pressure injector.

- Wear mask and sterile gloves.
- Practice a local disinfection of the implantation area (**with cleansing swabs containing 70% Isopropyl Alcohol, Dry**)
- Needle priming: with a Luer Lock syringe filled of the solution to inject (or sterile physiological serum), make a complete rinsing of the extension line in order to expel the air.
- Holding the needle from the sides, or the grip, remove the needle cap.
- Localize the port by palpation.
- Maintain firmly the port between 2 fingers.
- Prick between the fingers on centre of the septum from perpendicular to the skin (the septum is round)..
- Cross slowly the skin, then the septum (characteristic sensation as if it pricks a rubber) until the needle touch the port bottom (about 10 mm + peel thickness).
- Avoid the contact between the Huber needlepoint and the port bottom (risk of septum damage during the removal).**Remove the ONCO-GRIP® grip
- Fix the needle on the skin with a hypoallergenic dressing (preferably washable and transparent). **For the necessary checks after the introduction of the needle, see below.**
- Inject slowly and smoothly the solution or the physiological serum (rinsing).
- End of the procedure, rinse and lock the port according to the instructions of the manufacturer of the implantable port.
- Take off frankly the needle pushing the syringe piston.

VERIFICATIONS AFTER NEEDLE INTRODUCTION :

Remove the needle in one hand and simultaneously use the other hand to push the plunger of the syringe in order to apply a positive pressure.

WARNING: If you feel a resistance at the injection or if the solution doesn't inject it, don't force: this could be the sign of a port catheter thrombosis. Check firstly that the clamp is open and if the line is not bent.

Inform the responsible doctor and prepare a heparinisation of the port (according to the doctor's decision).

WARNING: Fully tighten all connections before use. Failure to attach caps can result in an embolism or bleeding.

WARNING: regularly check for leakage on the infusion line and connections.

WARNING: do not pull too strongly on the tubing before, during or after use. Too much pull strength applied may induce leakage.

INJECTION OF SEVERAL DIFFERENT SOLUTIONS :

WARNING: Never inject one after the other 2 products, which could precipitate or interact, without an intermediate rinsing with sterile physiological serum (thrombosis risk). In the event of doubt, make a rinsing.

WARNING: check for leaks

1. Close the clamp,

2. Take off the first syringe,

3. Adapt the second syringe,

WARNING: Take cares of not introduce air or infections during the change.

4. Open the clamp.

BLOOD SAMPLES :

1. Take 5 ml of blood and throw this sampling (because the blood is mixed to the heparin bolt)

2. Make the blood sampling.

3. Rinse immediately with 10 ml of physiological serum.

4. Heparinise

RINSING :

Injection of 3 to 5 ml of sterile physiological serum.

Directing the Huber point needle opening in the diametrically opposite direction of the implantable port exit channel increases the flushing efficiency.

HEPARINISATION :

Injection of 3 to 5 ml of a solution containing heparin (100 U/ml).

Routine heparinisation (follow the doctor's opinion)

- arterial catheter: every 7-15 days

- venous catheter: every 15-30 days

CONTRA-INDICATIONS :

On the skin against the septum:

• erythema, exanthema,

• unusual pain,

• weeping, discharge of fluid, fistula, abscess, any sign or suspected local infection,

• any sign or suspected turning around of the reservoir,

• difficulty or inability to locate the septum of the port.

Avoid direct contact between the device and alcohol-based solutions.

Solutions containing **Cabazitaxel** as active principle can not be injected with this needle.

Some solvents have been identified as being incompatible with a high concentration: **aromatic or halogenated hydrocarbons, esters, ethers and ketones.**

An adsorption phenomenon was observed on the following compounds:

- Bupivacaine hydrochloride	- Nicardipine hydrochloride
- Calcitriol	- Nifedipine
- Carmustine	- Nitroglycerin
- Ciclosporin	- Pentamidine isetonate
- Diazepam	- Propofol
- Fentanyl citrate	- Sufentanil citrate
- Heparin sodium	- Triflupromazine hydrochloride
- Insulin	- Voriconazole
- Isosorbide dinit	

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