



CATHÉTER D'EXTRACTION BILIAIRE BILIARY EXTRACTION CATHETER

3 LUMIÈRES

- 3 VRAIES LUMIÈRES POUR UNE PROCÉDURE SIMPLIFIÉE
- GRANDE RÉSISTANCE DU BALLONNET QUI PERMET UNE EXTRACTION DES CALCULS BILIAIRES EN TOUTE SÉCURITÉ
- BAGUE RADIO OPAQUE



- 3 LUMEN FOR A SIMPLIFIED PROCEDURE
- GREAT RESISTANCE OF THE BALLOON WHICH ALLOWS SAFETY EXTRACTION OF GALLSTONES
- RADIOPAQUE RING



Plusieurs variantes pour répondre à vos attentes :

- Le fil guide peut être positionné sur toute la longueur du cathéter (guidage long), ou émerger par un orifice situé à proximité de la partie distale (guidage distal). L'injection peut être faite avant ou après le ballonnet.
- 2 gammes de diamètres de gonflage possibles
- Egalement disponible en latex free

Several alternatives to answer your needs:

- The guide wire can be positioned over the entire length of the catheter (long guidance), or emerge by an opening located near the distal part (distal guidance). The injection can be made before or after the balloon.
- 2 ranges of inflation diameters
- Latex Free also available

3 LUMEN